

2016 Catalog

Oxford Instruments X-Ray Technology

Power Supplies, Integrated Sources and X-ray Tubes



Shielded X-ray
Tubes and
Power Supplies



Radiation Shielded
X-ray Tubes



Microfocus X-ray
Sources



Glass X-ray Tubes
with Window



Compact Integrated
X-ray Sources



Integrated X-ray
Sources

April / 2016

OXFORD
INSTRUMENTS

The Business of Science®

2016 Catalog

April 2016

Oxford Instruments X-Ray Technology is a global leader in the design and manufacture of integrated X-ray solutions. Our products range from integrated X-ray sources to high voltage power supplies and individual X-ray tubes featuring high stability, high X-ray flux and small spot sizes.

For more than 30 years, Oxford Instruments has been the best choice for analytical, industrial and medical original equipment manufacturers who demand the highest quality X-ray solutions.

This catalog should be used as a general guide for our products. The drawings in the datasheets are typical. Please ask us for the exact outline drawings and detailed specifications of the products.

If you do not see the product you are looking for in this catalog, please ask us about it!



ISO 9001:2008 Certified and RoHS Compliant



Solving Customer Specific Issues Since 1980

2016 Catalog

Table of Contents

Integrated X-ray Sources

Scafell Pike Integrated X-ray Source	4
Trinity 80kV, 33 μ m, 40 Watt Integrated X-ray Source	6
UltraBright 96000 Series 90kV Microfocus X-ray Source	8
Nova 96000 Series 90kV Water-Cooled Microfocus X-ray Source	12

Radiation Shielded X-ray Tubes

Pinnacles 50kV Microfocus X-ray Source	16
Neptune 5200 Series Water-Cooled Radiation Shielded X-ray Tube	18
Jupiter 5000 Series Radiation Shielded X-ray Tube	20
Apogee 5500 Series Radiation Shielded X-ray Tube	22

Potted X-ray Tubes

3000 Series 30kV X-ray Tube	24
-----------------------------	----

Glass X-ray Tubes

1000 Series Glass X-ray Tube	26
1500 Series Glass X-ray Tube	28
1501 Series Glass X-ray Tube	30
1550 Series Glass X-ray Tube	32

Power Supplies

Shasta Series X-ray Tube Power Supply	34
---------------------------------------	----

Appendix

Complete Listing of X-ray Tubes, Integrated Sources and Cables	36
Conditioning Procedure Application Note	42
X-ray Spectra Application Note	43
Operating Range Application Note	46
X-ray Tube Packaging Application Note	48
Beryllium Window Application Note	49
Filament Life Application Note	50
Shielding Application Note	52
Heat Management Application Note	55
X-ray Fluorescence Periodic Table	57



Scafell Pike is a compact and lightweight integrated X-ray source designed for handheld and portable XRF applications

Scafell Pike is a 50kV, 4W integrated X-ray source with a revolutionary stainless steel X-ray tube design leading to unprecedented ruggedness and flux stability in an incredibly light package. Scafell Pike is ideally suited to the unique challenges of handheld XRF and other applications where ruggedized design and close-coupled X-ray geometry are important.







Features	Benefits
Stable X-ray output	Highly repeatable measurements
Intelligent controller	Enables maximum filament life and minimizes startup times
Advanced digital interface	Simplifies communication and system integration
Small, compact design	Easily integrates into your handheld or benchtop X-ray systems
Fully shielded package	Eliminates X-ray leakage

Applications

- XRF
- Inspection
- Mining and geology
- Densitometry
- Thickness gauging
- Regulatory (RoHS/WEEE)
- Art and archeometry
- R&D

Specifications	
Operating Voltage Range:	4kV - 50kV
Maximum power:	4W
Beam current:	0-200µA
Focal spot size (maximum):	<800µm
Focus to Object Distance (FOD):	<1mm
Target material:	Rh
Window material and thickness:	Be, 127µm
Input voltage:	6.4 - 14 VDC
Input voltage ripple:	±5%
Control signals:	Digital communication, X-ray enable, and more (see next page)
Time to warm up to 50kV and 0.1mA:	<1 sec
Continuous flux stability:	<0.5%
Cooling method:	Conduction
Max case operating temp:	70°C
Ambient operating temp:	-10°C to 60°C
Radiation leakage:	<1.0µSv.hr ⁻¹ @ 10cm (@ 50kV, 80µ)
Weight:	<325g
Storage conditions:	Storage conditions:-25°C to 85°C* *Note: Humidity: 10-95% (no condensation) Condensation on Be window will cause window corrosion, vacuum loss, and X-ray tube failure

Scafell Pike Integrated X-ray Source



DIMENSIONS: [mm]
Inches

PIN	SIGNAL	VOLTAGE
1	SCL	0V - 5V
2	SDA	0V - 5V
3 - 4	INPUT POWER	6V - 18V
5 - 6	GND	0V
7 - 8	NO CONNECT	N/A
9	FILAMENT READY	0V - 5V
10	X-RAY ENABLE	0V - 5V
11	NO CONNECT	N/A

80kV, 33um, 40 Watt Integrated X-ray Source

Trinity

Trinity is an 80kV, 33µm, 40W, fully shielded integrated X-ray source designed for high resolution imaging applications, making it ideal for industrial inspection and non-destructive testing applications including PCB assembly, battery, plastic, metal and mechanical parts inspection, yet is versatile enough for medical imaging applications.

Trinity is configured in a compact, programmable package that simplifies communications and system integration without the need for high voltage cable connections or additional shielding, making it extremely reliable and cost-efficient. Side window and end window orientations are available.



Benefits

- Small, stable spot delivers distortion free measurements
- Compact, programmable, integrated package enables ease of installation and improved reliability
- XT proprietary shielding minimizes weight while ensuring extremely low radiation leakage
- LED indicators provide constant system status
- Fully lead-free and RoHS compliant


Applications

- Medical imaging
- Inspection of printed circuit boards and electronic devices
- Nondestructive testing of plastic, metal and mechanical parts
- CT imaging for life sciences and industrial inspection applications
- Food or Packaging imaging


Specifications	
Operating Voltage Range:	20-80kV
Maximum Power:	40W
Maximum Beam Current:	10-500µA
Target Material:	W
Focal Spot Size:	33µm (Nominal CEI/IEC 60336:2005)
Cone of Illumination:	37°
Spot to Window Spacing (FOD):	27mm ± 1mm
Window Material & Thickness:	1.40 mm (0.055") Glass & 3.60mm (0.142") Polystyrene
Accuracy:	kV: < 0.50% µA: < 1.0%
Flux Stability:	≤ 0.2% over 4-hour period
Rise Time:	≤ 200 ms (from standby)
Duty Cycle:	Continuous
Ripple:	1% RMS @ 80kV, 500µA
Temperature Coefficient:	1000ppm/°C
Temperature Conditions:	Operating: Maximum 55°C at tube heat sink Storage: -10°C to 50°C
Humidity:	15°C to 40°C ; 0-95% RH up to 5,000ft
Method of Cooling:	External cooling required and directed at the unit at 100 CFM
Thermal Cut-Off:	55°C ± 3°C at tube heat sink
Shielding:	Less than 0.1 mR/hour at 5cm away from surface of the unit when operated at 80kV, 500µA as per FDA 21 CFR 1020.40
Size:	10.5" L x 5.8" W x 4.0" H (266mm L x 146mm W x 101mm H)
Weight:	≤ 9 lbs.
Input Power:	22-26 VDC, 3A
Interface:	Analog (0-10VDC), RS422 & RS485
Safety & Regulatory Compliance:	Designed to meet UL, CE and RoHS Directive 2011/65/EU




Trinity 80kV, 33μm, 40W integrated X-ray source



The views above are consistent for both the end and side window configurations.



End Window Configuration Part Number 9100005



Side Window Configuration Part Number 9100004

DIMENSIONS: Inches [mm]

Caution:
*M5 mounting holes have a maximum screw penetration of 6mm.
*M6 mounting holes have a maximum screw penetration of 8mm.

Control Interface Box PN 9100006: An optional adapter for manual analog control may be ordered – useful for lab use and qualification testing

90kV Microfocus X-ray Source

UltraBright 96000 Series



The UltraBright Microfocus System 96000 Series is a 90kV, 80W X-ray source designed for applications where high brightness, high magnification and small spot size are important.

Operated by an external high voltage Smart Controller capable of providing variable voltage and power control, the UltraBright Microfocus System delivers exceptional magnification and image quality with full control of “Brightness”. Maximum flux output is maintained through automatic matching of a given power setting to a corresponding optimal spot size.



Benefits

- Exceptional magnification and image quality
- High power operation — ideal for high flux applications and experiments.
- Integrated package eliminates HV cable for improved reliability
- Complete range of user control — ideal for research applications
- Compact, lightweight design — ideal for portable applications

Applications

- Microtomography
- Microdiffraction
- Microfluorescence
- CT imaging for life sciences and industrial inspection


Product Ordering Table

Target Material (Part#)	Voltage	Power	Power Density
W (96004)	10-90kV	10-80W	2.5W/μm
Mo (96002)	20-60kV	20-60W	1.5 W/μm
Cu (96000)	20-60kV	20-60W	1.5 W/μm


Specifications	
Operating voltage range:	See product ordering table
Maximum Power:	See product ordering table
Maximum beam current:	2.0mA
Focal spot size:	14-20μm @ maximum voltage and minimum power
Focus to Object Distance (FOD):	4mm
Cone of illumination:	50° x 74° (nominal) See chart on next page
Window material and thickness:	Be, 254μm
Window diameter (unobstructed):	9.5mm (0.37")
Window configuration	End window
Target material:	See product ordering table
Ambient operating temperature:	10°C to 40°C
Maximum operating temperature (anode):	70°C
Cooling method:	Forced air (150 CFM @ 4" recommended for continuous operation)
Shielding:	Not shielded
Dimensions:	392.4mm L x Ø114.3mm (15.5" L x Ø4.5")
Weight:	≤ 4kg (8.81lbs.)
Storage conditions:	-10°C to 55°C Barometric Pressure: 50-106kPa; Humidity: 10-90% (no condensation) Condensation on Be window will cause window corrosion, vacuum loss and X-ray tube failure.




UltraBright 96000 Series 90kV Microfocus X-ray Source



DIMENSIONS: [mm]
Inches



Radiation Pattern as Seen Normal to Window



Dimension	Description	Units	Emitted Cone and Spot Position		
			Farthest	Nominal	Nearest
A	Location of radiation cone center	Degrees	10.6	12.9	16.4
B	Radius of cone	Degrees	32.2	36.8	42
C	Window to spot distance	mm	4.47	3.14	1.82



The source (left) is connected to the controller (right) with a DB-25 cable.

The female end of the DB-25 cable connects to the source and the male end of the DB-25 cable connects to the controller.



90kv Microfocus X-ray Source

UltraBright 96000 Series (cont.)



RS232 Control Command Set

Protocol: RS-232-C
 Baud Rate: 9600 ASYNC
 Flow control: None
 Data bits: 8
 Stop bits: 1
 Parity: None
 Connector: Type: 25 pin

Controller Unit Specifications	
Functions:	Key switched power, HV on/off, kV adjust, brightness/autofocus adjust
External control:	Remote control
Power consumption:	100W maximum
Input voltage:	110/240 AC autosensing
Approximate weight:	4kg
HV cable:	Not necessary
LV cable:	Std 25 pin D-type connector (15 feet long)

Functions


Anode voltage set: 10 to 90kV (example: VCN 50 = set 50kV)
 Brightness set: 10 to 80W (example: WCN 40 = set 40W)
 Command: X-ray ON/OFF
 Command: Voltage min-max set
 Command: Brightness min-max set

Item	RS232 Control	Manual Operation	Notes
Remote/local switch	No	Yes	Switch is located on rear panel
Power on/off	No	Yes	For remote operation, front panel on/off switch must be "On"
X-ray on/off:	Yes	Yes	For remote operation, front panel on/off switch must be "On"
Voltage up/down	Yes	No	Front panel switch disabled
Brightness control	Yes	No	Front panel switch disabled

Read Back


Voltage: (example: VM 30 = 30kV)
 Brightness: (example: WM 20 = 20W)
 Status: Stand-by, warm-up, output, fault modes in ASCII format
 Fault: Display panel information except remote/local mode will be in ASCII format
 Other: ROM version number

Power on can be accomplished remotely by X-ray on/off command. However, if cathode emitter is turned off, power is restored only through front panel on/off switch. Safety interlocks available on real panel.



Controller Unit Reference Drawing / Dimensions in mm

Software Control Option



- Runs under LabVIEW RT & MS-Windows
- Works in conjunction with RS232 control interface
- Complete control of voltage, power, and focus
- Dynamic status display
- Dynamic fault display
- Ideal for R&D applications
- Open software architecture allows for modification with available additional development software



Functions	Software Control	Manual Operation	Notes
Remote/Local switch	No	Yes	Switch is located on rear panel
Power On/Off	No	Yes	For software operation, front panel On/Off switch must be "On"
X-ray On/Off	Yes	Yes	For software operation, front panel On/Off switch must be "On"
Voltage up/down	Yes	No	Front panel switch disabled
Power adjust	Yes	No	Front panel Brightness dial disabled

Microfocus Source with the following Class-Leading Performance Characteristics:

- The UltraBright 96000 Series is a fully integrated 90kV X-ray source. Its high voltage power supply and controller provide variable control of high voltage from 10-90kV and beam current from .33-2 milliamps with full control of "Brightness". The Smart Controller calculates spot size for a given power setting for maximum flux output.
- Voltage and current rating (90kV, 2.0mA) are subject to maximum power dissipation rating of 80W. The X-ray tube assembly is sealed, air-cooled, and rated for continuous operation.
- X-ray microfocus spot size is continuously adjustable from 14 μ m to 20 μ m. Power de-rating is provided at small spot sizes but source power is greater than or equal to 20W for a 20 μ m spot size.
- The anode target material is comprised of Tungsten as standard, however other targets are available (Cu, Mo). The target is inclined at a takeoff angle of 15 degrees with respect to the electron beam, and the exit window is aligned at an angle of 30 degrees with respect to the electron beam, so that a round microfocus X-ray spot is projected through the exit window.
- The stability of the microfocus X-ray spot shall be less than 5 μ m RMS over a period of 8 hours, as verified by test. A warm-up time of up to two hours is necessary in order to meet this specification.
- The system is supplied with a 254 micron Be exit window, allowing for close coupling (4mm) of object with the anode X-ray spot.
- LabVIEW RT Software Interface: The Smart Controller is outfitted with a software package that provides remote control of the various functions, such as kV, mA, Brightness, power etc. It includes an RS232 Communication package and an RT version of National Instruments LabVIEW. See Software control datasheet for complete description.

90kV Water-Cooled Microfocus X-ray Source

Nova 96000 Series



The Nova Microfocus System 96000 Series is a 90kV, 80W, water-cooled X-ray source designed for applications where high power, high magnification and small spot size are important.

Operated by an external high voltage Smart Controller capable of providing variable voltage and power control, the Nova Microfocus System delivers exceptional magnification and image quality with full control of “Brightness.” Maximum flux output is maintained through automatic matching of a given power setting to a corresponding optimal spot size.



Benefits

- Exceptional magnification and image quality
- High power operation — ideal for high flux applications and experiments.
- Integrated package eliminates HV cable for improved reliability
- Complete range of user control — ideal for research applications
- Compact, lightweight design — ideal for portable applications

Applications

- Microtomography
- Microdiffraction
- Microfluorescence
- CT imaging for life sciences and industrial inspection


Product Ordering Table

Target Material (Part#)	Voltage	Power	Power Density
W (96013)	10-90kV	10-80W	2.5W/μm
Mo (96016)	20-60kV	20-60W	1.5 W/μm


Specifications	
Operating voltage range:	See product ordering table
Maximum power:	See product ordering table
Maximum beam current:	2.0mA
Focal spot size:	14-20μm @ maximum voltage and minimum power
Focus to Object Distance (FOD):	4mm
Cone of illumination:	50° x 74° (nominal) See chart on next page
Window material and thickness:	Be, 254μm
Window diameter (unobstructed):	9.5mm (0.37")
Window configuration	End window
Target material:	See product ordering table
Ambient operating temperature:	10°C to 40°C
Maximum operating temperature (anode):	70°C
Cooling method:	H ₂ O 0.15 l/min @ 15 psi
Shielding:	Not shielded
Dimensions:	392.4mm L x 114.3mm W (15.5" L x 4.5" W)
Weight:	≤ 4kg (8.81 lbs)
Storage conditions:	-10°C to 55°C Barometric Pressure: 50-106kPa; Humidity: 10-90% (no condensation) Condensation on Be window will cause window corrosion, vacuum loss and X-ray tube failure



Nova 96000 Series 90kv Water-Cooled Microfocus X-ray Source



Dimension	Description	Units	Emitted Cone and Spot Position		
			Farthest	Nominal	Nearest
A	Location of radiation cone center	Degrees	10.6	12.9	16.4
B	Radius of cone	Degrees	32.2	36.8	42
C	Window to spot distance	mm	4.47	3.14	1.82



The source (left) is connected to the controller (right) with a DB-25 cable.



The female end of the DB-25 cable connects to the source and the male end of the DB-25 cable connects to the controller.

90kv Water-Cooled Microfocus X-ray Source Nova 96000 Series (cont.)

RS232 Control Command Set

Protocol: RS-232-C
 Baud Rate: 9600 ASYNC
 Flow control: None
 Data bits: 8
 Stop bits: 1
 Parity: None
 Connector: Type: 25 pin

Functions

Anode voltage set: 10 to 90kV (example: VCN 50 = set 50kV)
 Brightness set: 10 to 80W (example: WCN 40 = set 40W)
 Command: X-ray ON/OFF
 Command: Voltage min-max set
 Command: Brightness min-max set


Read Back

Voltage: (example: VM 30 = 30kV)
 Brightness: (example: WM 20 = 20W)
 Status: Stand-by, warm-up, output, fault modes in ASCII format
 Fault: Display panel information except remote/local mode will be in ASCII format
 Other: ROM version number

Controller Unit Specifications	
Functions:	Key switched power, HV on/off, kV adjust, brightness/autofocus adjust
External control:	Remote control
Power consumption:	100W maximum
Input voltage:	110/240 AC autosensing
Approximate weight:	4kg
HV cable:	Not necessary
LV cable:	Std 25 pin D-type connector (15 feet long)


Item	RS232 Control	Manual Operation	Notes
Remote/local switch	No	Yes	Switch is located on rear panel
Power on/off	No	Yes	For remote operation, front panel on/off switch must be "On"
X-ray on/off:	Yes	Yes	For remote operation, front panel on/off switch must be "On"
Voltage up/down	Yes	No	Front panel switch disabled
Brightness control	Yes	No	Front panel switch disabled

Power on can be accomplished remotely by X-ray on/off command. However, if cathode emitter is turned off, power is restored only through front panel on/off switch. Safety interlocks available on real panel.



Controller Unit Reference Drawing / Dimensions in mm

Software Control Option



- Runs under LabVIEW RT & MS-Windows
- Works in conjunction with RS232 control interface
- Complete control of voltage, power, and focus
- Dynamic status display
- Dynamic fault display
- Ideal for R&D applications
- Open software architecture allows for modification to user interface with available additional development software



Functions	Software Control	Manual Operation	Notes
Remote/Local switch	No	Yes	Switch is located on rear panel
Power On/Off	No	Yes	For software operation, front panel On/Off switch must be "On"
X-ray On/Off	Yes	Yes	For software operation, front panel On/Off switch must be "On"
Voltage up/down	Yes	No	Front panel switch disabled
Power adjust	Yes	No	Front panel Brightness dial disabled

Microfocus Source with the following Class-Leading Performance Characteristics:

- The Nova 96000 Series is a fully integrated 90kV X-ray source. Its high voltage power supply and controller provide variable control of high voltage from 10-90kV and beam current from .33-2 milliamps with full control of "Brightness". The Smart Controller calculates spot size for a given power setting for maximum flux output.
- Voltage and current rating (90kV, 2.0mA) are subject to maximum power dissipation rating of 80W. The X-ray tube assembly is sealed, water-cooled, and rated for continuous operation.
- X-ray microfocus spot size is continuously adjustable from 14 μ m to 20 μ m. Power de-rating is provided at small spot sizes but source power is greater than or equal to 20W for a 20 μ m spot size.
- The anode target material is comprised of Tungsten as standard, however a molybdenum target is also available. The target is inclined at a takeoff angle of 15 degrees with respect to the electron beam, and the exit window is aligned at an angle of 30 degrees with respect to the electron beam, so that a round microfocus X-ray spot is projected through the exit window.
- The stability of the microfocus X-ray spot shall be less than 5 μ m RMS over a period of 8 hours, as verified by test. A warm-up time of up to two hours is necessary in order to meet this specification.
- The system is supplied with a 254 micron Be exit window, allowing for close coupling (4mm) of object with the anode X-ray spot.
- LabVIEW RT Software Interface: The Smart Controller is outfitted with a software package that provides remote control of the various functions, such as kV, mA, Brightness, power etc. It includes an RS232 Communication package and an RT version of National Instruments LabVIEW. See Software control datasheet for complete description.

50kV Microfocus X-ray Source Pinnacles 50kV

Developed for applications that require high resolution over a wide-angle field of view, the Pinnacles 50kV Microfocus X-ray source features high flux output.

Its compact design is fully radiation shielded and insulated with an integrated high voltage cable located on the side of the tube for easy connection.

The Shasta μ F power supply has been optimized to power the Pinnacles 50kV Microfocus X-ray tube.



Benefits

- Wide operating range enables optimal image contrast
- Wide field of view
- Fully shielded package eliminates X-ray leakage and easily integrates into your system
- Integrated high voltage cable

Applications

- Medical imaging
- Printed circuit board and electronic device inspection
- Nondestructive testing of plastic, metal and mechanical parts

Specifications	
Operating voltage range:	10-50kV
Maximum power:	12W
Maximum beam current:	1.0mA
Focal spot size:	10 μ m (50kV, 3W) line pair resolution using JIMA RT RC-02
Focus to Object Distance (FOD):	35.18mm (1.385")
Target material:	W
Target angle:	45°
Cone of illumination (unobstructed):	40.5° \pm 0.5°
Window material and thickness:	Be, 254 μ m
Window diameter (unobstructed):	16.88mm (.66")
Maximum operating temperature:	50°C at potting surface
Ambient operating temperature:	0°C to 40°C; 0-95% RH up to 5,000ft
Cooling method:	Forced air @ 150cfm at 100mm (4.0") recommended
Shielding:	Fully shielded. X-ray leakage < 1.0 μ Sv.hr-1 at 10cm
Weight:	\approx 1.37kg (3 lbs)
Storage conditions:	-10°C to 55°C; Barometric Pressure: 50-106kPa; Humidity: 10-90% (no condensation) Condensation on Be window will cause window corrosion, vacuum loss and X-ray tube failure


Shasta μ F Power Supply

- Industry-standard 24V Input
- High voltage, cathode, and grid controls
- Intuitive analog control interface
- Focusing grid adjustment for optimum spot size




Pinnacles 50kV Microfocus X-ray Source

Pinnacles Ordering Table	
Part Number	High Voltage Cable Length
9400001	18 inch
9400002	39 inch [1 meter]
9400003	79 inch [2 meter]
9400004	118 inch [3 meter]



Low Voltage Cable Ordering Table	
Part Number	Cable Length
9290002	39 inch [1 meter]
9290003	79 inch [2 meter]
9290004	118 inch [3 meter]

LV Connector Pin	Tube Pin
3	GROUND
4	HEATER
5	HEATER GROUND
7	L1
12	L2
14	L3



Notes

1. THIS TUBE IS FULLY RADIATION SHIELDED TO 50kV/12W EXCEPT 40° X-RAY CONE.
2. DIMENSIONS ARE IN INCHES. DIMENSIONS [] ARE IN MILLIMETERS.

Water-Cooled Radiation Shielded X-ray Tube

Neptune 5200 Series

X-RAY

The Neptune 5200 Series is a water-cooled 50kV, 100W packaged X-ray tube designed for applications where high flux density and continuous operation are important.

Utilizing our high stability and high intensity X-ray tube technology, the Neptune 5200 Series is ideal for most industrial inspection and non-destructive testing applications that require high resolution, including plastic, metal and mechanical parts inspection. Flexible and reliable, this unit is also highly suited for use in high power XRF applications.

The 5200 Series has a brass package that utilizes 0.2 liter/min of water flow, which enables the unit to provide maximum X-ray shielding and heat dissipation. The design includes high voltage, filament and water flow connectors, making it ideal for plug and play operation.

The Neptune 5200 Series is available in wide range of targets and price points to meet your needs.



Benefits

- Wide operating range enables optimal image contrast
- Stable X-ray output delivers high precision measurements
- Low attenuation beryllium window ensures high transmission of low energy X-rays
- Fully-shielded compact package eliminates X-ray leakage and easily integrates into your system


Applications

- Non-destructive testing of plastic, metal and mechanical parts
- Thickness gauging
- Analytical XRF

Specifications	
Operating Voltage Range:	10-50kV
Maximum Power:	100W
Maximum Beam Current:	2.0mA
Maximum Filament Current:	2.40A
Filament Voltage:	3.75V (Nominal)
Target Material:	See Product Ordering Table
Spot Size:	175 μ m where $X+Y/2$ and $X < 210\mu$ m and $Y < 210\mu$ m
Cone of Illumination:	25°
Spot to Window Spacing (FOD):	48.8 mm \pm 1mm (1.92")
Window Material & Thickness:	Be @ 127 μ m
Flux & Current Stability:	\leq 0.2% over 4-hour period
Duty Cycle:	Continuous
Ambient Temperature Conditions:	Operating: 0 to 40°C Storage: -10°C to 50°C
Humidity:	0-95% RH up to 5,000ft
Method of Cooling:	Water cooling $>$.21 l/min. Forced air cooling directed at the unit at 150 CFM may be required at high power operation. Must not exceed 55°C at case surface.
Shielding:	0.25mR/hr @ 2" (except HV connection through HV cable)
Dimensions:	210mm L X 106 mm W (8.25" L X 4.18" W)
Weight:	6.17 kg (13.6 lbs)



Neptune 5200 Series Water-Cooled Radiation Shielded X-ray Tube



Product Ordering Table	
Part Number	Target Material
93221	Rh
93211	Mo
93212	W

The Jupiter 5000 Series is a 50kV, 50W packaged X-ray tube designed for applications where high flux density and continuous operation are important.

Utilizing our highly stable and high intensity X-ray tube technology, the Jupiter 5000 Series is ideal for medical imaging applications and most industrial inspection and non-destructive testing applications that require high resolution, including PCB assembly, battery, plastic, metal and mechanical parts inspection.

The 5000 Series features a stainless steel, lead-lined package that is filled with dielectric oil, which enables the unit to provide maximum X-ray shielding and heat dissipation. The design includes high voltage and filament connectors, making it ideal for plug and play operation.

The Jupiter 5000 Series is available in a wide range of spot sizes, targets and price points to meet your needs.



Benefits

- Wide operating range enables optimal image contrast
- Stable X-ray output delivers high precision measurements
- Low attenuation beryllium window ensures high transmission of low energy X-rays
- Fully-shielded compact package eliminates X-ray leakage and easily integrates into your system


Applications

- Medical Imaging
- Printed circuit board and electronic device inspection
- Non-destructive testing of plastic, metal and mechanical parts
- Thickness gauging
- Analytical XRF

Specifications

Operating Voltage Range:	10-50kV
Maximum Power:	50W
Maximum Beam Current:	1.0mA
Focal spot size:	P/N 93089: 50µm P/N 93095: 55µm
Maximum Filament Current:	1.7A
Filament Voltage:	2.0V (nominal)
Focus to Object Distance (FOD):	See diagram next page
Window material and thickness:	Be @ 127µm
Cone of illumination (unobstructed):	23°
Window diameter (unobstructed):	11.43mm (.450")
Target material:	See product ordering table next page
Target angle:	12°
Stability:	0.2% 4 hours
Polarity:	Grounded cathode
Maximum operating temperature:	55°C on case surface
Ambient operating temperature:	0°C to 40°C
Cooling method	Forced air @ 150cfm
Shielding:	0.25mR/hr @ 2" (except at HV connection)
Dimensions:	180mm L x Ø70mm (7.09" L x Ø2.76")
Weight:	2.26kg (5.0 lbs)
Storage Conditions:	-10°C to 55°C Barometric Pressure: 50-106kPa; Humidity: 10-90% (no condensation) Condensation on Be window will cause window corrosion, vacuum loss and X-ray tube failure

Jupiter 5000 Series Radiation Shielded X-ray Tube



Notes

- 1. USE AN AS568A-139 O-RING IN MOUNTING PLATE GROOVE IF REQUIRED
- 2. DIMENSIONS ARE IN INCHES. DIMENSIONS [] ARE IN MM

Product Ordering Table		
Target Material	Voltage	Power
W (93089)	10-50kV	50W
Mo (93095)	10-50kV	50W



Use a Shasta Power Supply for Peak Performance

The Apogee 5500 Series is a 50kV, 50W packaged X-ray tube designed for applications where high flux density and continuous operation are important.

Utilizing our high stability, high intensity X-ray tube technology coupled with grid-controlled variable focus enables our Apogee design to produce very small focal spots; this makes the Apogee 5500 Series ideal for most industrial inspection and non-destructive testing applications that require high resolution, including PCB assembly, battery, plastic, metal and mechanical parts inspection. Flexible and reliable, this unit is also well suited for use with X-ray optics.

The Apogee 5500 Series is configured in a compact stainless steel, lead-lined package filled with dielectric oil, which enables the unit to provide maximum X-ray shielding and heat dissipation. The design includes high voltage and filament connectors, making it ideal for plug and play operation.



Benefits


- Wide operating range enables optimal image contrast
- Stable X-ray output delivers high-precision measurements
- Low attenuation beryllium window ensures high transmission of low energy X-rays
- Fully shielded compact package eliminates X-ray leakage and easily integrates into your system

Applications

- Medical Imaging
- Inspection of printed circuit boards and electronic devices
- Non-destructive testing of plastic, metal and mechanical parts
- Thickness gauging
- Analytical XRF

Specifications	
Operating Voltage Range:	10-50kV
Maximum Power:	50W
Maximum Beam Current:	1.0mA
Grid Voltage:	0-100V (Oxford Shasta Power Supply recommended)
Maximum Filament Current:	1.7A
Filament Voltage:	2.0V (Nominal)
Target Material:	See product ordering table on next page
Focal Spot Size:	35µm *nominal per IEC60336, NEMA XR5-1992 (R1999)
Cone of Illumination:	22°
Spot to Window Spacing (FOD):	31.05mm ±1mm
Window Material and Thickness:	Be @ 127µm
Flux and Current Stability:	≤ 0.2% over 4-hour period
Duty Cycle:	Continuous (150 CFM airflow required)
Ambient Temperature Conditions:	Operating: 0 to 40°C Storage: -10°C to 50°C
Humidity:	0-95% RH up to 5,000 feet
Cooling Method:	External cooling required and directed at the unit at 150 CFM. Must not exceed 55°C at case surface.
Shielding:	0.25mR/hr @ 2" (except HV connection through HV cable)
Dimensions:	180mm L X Ø70mm (7.09" L X Ø2.76")
Weight:	1.82kg (4.0 lbs)

Apogee 5500 Series Radiation Shielded X-ray Tube



Pin	Oxford Cable Wire	Description
A	Red	Filament
B	Black	Filament Return
C	Clear	Grid Bias
Backshell	Eyelet	Ground

Product Ordering Table	
Part Number	Target Material
93501	W
93502	Mo
93504	Rh
93510	Cu



Use a Shasta
Power Supply for
Peak Performance

The Oxford Instruments 3000 Series X-ray tube has been developed for high flux stability and long life, making it ideal for continuous operation.

A low cost answer for high spectral purity radiation, the 3000 Series is encapsulated in silicone rubber and features a grounded cathode and low attenuation Beryllium window.

Features	Benefits
Continuous operation: 9W	High sensitivity and high precision measurement
Beryllium window	Higher flux of low-energy X-rays, especially from target L series lines
Compact, insulated light-weight package	Configuration allows flexible installation

Ordering Information	
Part Number	Target Material
90151	Rh
90152	Ag
90153	Mo




Applications

- Analytical (XRF)
- Particle Analysis
- Thickness Gauging
- Soft X-ray Radiography
- Spectroscopy
- Stress Analysis

Specifications	
Operating Voltage Range:	4-30kV
Maximum power:	9W
Maximum beam current:	0.3mA
Maximum filament current:	2.0A
Filament voltage:	1.75V (nominal)
Focal spot size:	1.0mm (nominal)
Focus to Object Distance (FOD):	28.2mm (1.1")
Target material:	See ordering table above
Window material and thickness:	Be @ 127µm
Unobstructed cone of illumination:	28°
Unobstructed window diameter:	10.4mm (.41")
Target angle:	20°
Shielding:	Partially radiation shielded (see drawing next page)
Weight:	260g
Cooling method:	Forced air: 150CFM @ 100mm (4.0") and appropriate heat sink recommended for full power
Maximum operating temp:	50°C at potting surface
Ambient operating temp:	0°C to 40°C
Storage conditions:	-40°C to 70°C*
	*Note: Barometric Pressure: 50-106kPa Humidity: 10-90% (no condensation) Condensation on Be window will cause window corrosion, vacuum loss, and X-ray tube failure

3000 Series 30kV X-ray Tube



Notes

1. WARNING: THIS TUBE IS NOT RADIATION SHIELDED.
2. DIMENSIONS ARE IN INCHES. DIMENSIONS [] ARE IN MM.

Oxford Instruments glass X-ray tubes are recognized for their performance and long life.

High flux and spot size stability make our X-ray tubes an ideal solution for demanding applications, such as those requiring continuous operation. The 90507 is uniquely designed with a very small isostatically focused spot for high resolution applications, such as mini C-Arm fluoroscopy. The robust electron gun assembly has been constructed for optimal use in integrated X-ray sources, where heat dissipation is an issue. Long tube life is achieved by ultra-high vacuum maintained with the Oxford Instruments unique Pin Flash getter. This tube operates in bi-polar mode.



Benefits


- Exceptional image quality
- Stable X-ray output delivers high precision measurements
- Small, stable spot delivers distortion-free measurements
- RoHS compliant design

Applications

- CT imaging for life sciences and industrial inspection
- Densitometry
- Thickness gauging
- Phase contrast imaging
- Medical imaging

Specifications	
Operating voltage range:	40-80kV (bi-polar operation: -40kV cathode, +40kV anode)
Maximum Power:	40W continuous
Maximum beam current:	0.5mA
Focal spot size:	33µm *nominal per IEC60336, NEMA XR5-1992 (R1999)
Maximum filament current:	1.7A
Filament voltage:	2.0V (nominal)
Focus to Object Distance (FOD):	14.2mm (0.56") (nominal)
Window material and thickness:	Glass—1.40mm ± 0.15
Target material:	W
Target angle:	20°
Maximum oil temperature:	80°C
Cooling method:	Oil
Weight:	114g (0.25lbs)
Storage conditions:	-10°C to 55°C Barometric Pressure: 50-106kPa; Humidity: 10-90% (no condensation)

1000 Series - 90507 Glass X-ray Tube



Notes

1. Dimensions are in inches.
2. This X-ray tube is designed to operate in an oil filled high voltage enclosure. Do not allow the oil to exceed 80°C . Proper operation of the X-ray tube requires cooling oil to circulate freely around the X-ray tube envelope.
3. This X-ray tube produces X-rays in all directions. As such, it must only be operated in a radiation-shielded enclosure.
4. Tubes to be shipped with two teflon-coated copper wire leads, 18 AWG X 12.0 MIN, soldered to pins #6 and #2.

The 1500 Series X-ray tube is a 50kV, 50W X-ray tube designed for applications where high flux density and continuous operation are important.

Utilizing our highly stable and high intensity X-ray tube technology, the 1500 Series X-ray tube is ideal for medical imaging, XRF applications and most industrial inspection and non-destructive testing applications that require high resolution, including PCB assembly, battery, plastic, metal and mechanical parts inspection.



The 1500 Series X-ray tube can also be supplied in a stainless steel, lead-lined package that is filled with dielectric oil that enables the unit to provide maximum X-ray shielding and heat dissipation, effectively replicating our popular Jupiter 5000 Series packaged tube, which includes high voltage and filament connectors making it ideal for plug and play operation.

The 1500 Series X-ray tube is available in a wide range of spot sizes, targets and price points to meet your needs.

Benefits


- Wide operating range enables optimal image contrast
- Stable X-ray output delivers high precision measurements
- Low attenuation beryllium window ensures high transmission of low energy X-rays

Applications

- Medical imaging
- Inspection of printed circuit boards and electronic devices
- Nondestructive testing of plastic, metal and mechanical parts
- Thickness gauging
- Analytical XRF

Specifications	
Operating Voltage Range:	10-50kV
Maximum Power:	50W
Maximum Beam Current:	1.0mA
Maximum Filament Current:	1.70A
Filament Voltage:	2.0V (Nominal)
Target Material:	All
Spot Size:	50-150µm (nominal per IEC60336,NEMA XR5-1999)
Cone of Illumination:	23°
Spot to Window Spacing (FOD):	30.8 mm ± 1mm (1.2")
Window Material and Thickness:	Be @ 127µm
Flux & Current Stability:	≤ 0.2% over 4-hour period
Duty Cycle:	Continuous
Ambient Temperature Conditions:	Operating: 0°C to 40°C Storage: -10°C to 50°C
Humidity:	0-95% RH up to 5,000ft
Method of Cooling:	Must not exceed 80°C oil temperature. Customer provides enclosure and cooling.
X-ray Shielding:	Customer must provide enclosure with adequate shielding. Tube emits X-rays in all directions.
Dimensions:	81mm L X 47mm W (3.2" L X 1.8" W)
Weight:	119g

1500 Series Glass X-ray Tube



DIMENSIONS: [mm]
Inches



Use a Shasta Power Supply for Peak Performance

Glass X-ray Tube 1501 Series

The 1501 Series X-ray tube is a 50kV, 50-75W X-ray tube designed for applications where high current, high flux density and continuous operation are important.

Utilizing our highly stable and high intensity X-ray tube technology, the 1501 Series X-ray tube is ideal for medical imaging, XRF applications and most industrial inspection and non-destructive testing applications that require high resolution, including PCB assembly, battery, plastic, metal and mechanical parts inspection.



The 1501 Series X-ray tube can also be supplied in a stainless steel, lead-lined package that is filled with dielectric oil that enables the unit to provide maximum X-ray shielding and heat dissipation.

The 1501 Series was designed in response to the need for higher current coupled with lower operating potentials.

The 1501 Series X-ray tube is available in a wide range of spot sizes, targets and price points to meet your needs.

Benefits

- Wide operating range enables optimal image contrast
- Stable X-ray output delivers high precision measurements
- Low attenuation beryllium window ensures high transmission of low energy X-rays


Applications

- Medical imaging
- Inspection of printed circuit boards and electronic devices
- Nondestructive testing of plastic, metal and mechanical parts
- Thickness gauging
- Analytical XRF

Specifications	
Operating Voltage Range:	4-50kV
Maximum Power:	50-75W
Maximum Beam Current:	2.5mA
Maximum Filament Current:	2.4A
Filament Voltage:	3.75V (Nominal)
Target Material:	Rh, Cr, Mo
Spot Size:	125µm (nominal per IEC60336, NEMA XR5-1999)
Cone of Illumination:	22° Minimum
Spot to Window Spacing (FOD):	30.8 mm ± 1mm (1.2")
Window Material and Thickness:	Be @ 127µm
Flux & Current Stability:	≤ 0.2% over 4-hour period
Duty Cycle:	Continuous
Ambient Temperature Conditions:	Operating: 0°C to 40°C Storage: -10°C to 50°C
Humidity:	0-95% RH up to 5,000ft
Method of Cooling:	Must not exceed 80°C oil temperature. Customer provides enclosure and cooling.
X-ray Shielding:	Customer must provide enclosure with adequate shielding. Tube emits X-rays in all directions.
Dimensions:	81mm L X 47mm W (3.2" L X 1.8" W)
Weight:	119g



1501 Series Glass X-ray Tube



DIMENSIONS: [mm]
Inches

The 1550 Series X-ray tube is a 50kV, 50W X-ray tube designed for applications where high flux density and continuous operation are important.

Utilizing our highly stable, high intensity X-ray tube technology coupled with grid-controlled variable focus enables our 1550 Series X-ray tube to produce very small focal spots; this makes the 1550 Series ideal for most industrial inspection and non-destructive testing applications that require high resolution, including PCB assembly, battery, plastic, metal and mechanical parts inspection. Flexible and reliable, this unit is also highly suited for use with X-ray optics.



The 1550 Series X-ray tube can also be supplied in a stainless steel, lead-lined package that is filled with dielectric oil that enables the unit to provide maximum X-ray shielding and heat dissipation; this configuration is our popular Apogee 5500 Series packaged tube, which includes high voltage and filament connectors making it ideal for plug and play operation.

Benefits



- Wide operating range enables optimal image contrast
- Stable X-ray output delivers high precision measurements
- Low attenuation beryllium window ensures high transmission of low energy X-rays

Applications

- Medical imaging
- Inspection of printed circuit boards and electronic devices
- Nondestructive testing of plastic, metal and mechanical parts
- Thickness gauging
- Analytical XRF

Specifications	
Operating Voltage Range:	10-50kV
Maximum Power:	50W
Maximum Beam Current:	1.0mA
Grid Voltage:	0-100V
Maximum Filament Current:	1.70A
Filament Voltage:	2.0V (Nominal)
Target Material:	Cu, W, Mo, Co, Rh
Spot Size:	<50µm (X and Y)
Cone of Illumination:	22°
Spot to Window Spacing (FOD):	30.8 mm ± 1mm (1.213")
Window Material and Thickness:	Be @ 127µm
Flux & Current Stability:	≤ 0.2% over 4-hour period
Duty Cycle:	Continuous
Ambient Temperature Conditions:	Operating: 0°C to 40°C Storage: -10°C to 50°C
Humidity:	0-95% RH up to 5,000ft
Method of Cooling:	Must not exceed 80°C oil temperature. Customer provides enclosure and cooling.
X-ray Shielding:	Customer must provide enclosure with adequate shielding. Tube emits X-rays in all directions.
Dimensions:	81mm L X 47mm W (3.2" L X 1.8" W)
Weight:	119g

1550 Series Glass X-ray Tube



DIMENSIONS: Inches



Use a Shasta Power Supply for Peak Performance

Oxford Instruments Shasta series power supply features a robust design that has been optimized to power grounded filament X-ray tubes from Oxford Instruments, yet its versatility enables it to power virtually any grounded filament X-ray tube.

Utilizing closed loop emission control circuitry that delivers low ripple, Shasta provides highly regulated beam current and high stability resulting in superior performance. Local and remote analog control enables convenient operation in setting voltage & emission current.

- Models with grid focus control are designed to provide optimal performance with our Apogee tubes

Benefits

- Compact Design
- Adjustable Emission Current
- Voltage & Current Programming
- Safety Interlock
- Bias Voltage Option Available
- UL, CE & TUV Certified

Applications



- XRF, XRD, Medical Imaging, Industrial Inspection & NDT



Specifications	
Operating voltage range:	0-50kV or 0-60kV models (see product ordering table next page)
Maximum Power:	50W or 60W models (see product ordering table next page)
Maximum beam current:	1.0mA
DC Filament Supply:	Current: 0.3 to 3.5A Voltage: 0 to 5.0 VDC
Voltage Regulation:	Load: < 0.01 % for 50% of max load variation Line: < 0.01% for 10% change in input voltage
Current Regulation:	Load: < ± 2µA (Beam Current) Line: < ± 2µA (Beam Current)
Ripple:	< 100V peak to peak
Stability:	± 0.1% over an 8-hour period after 30-minutes warm-up
Input Voltage & Power:	24VDC, ± 10%; 100 Watts
Voltage Control:	Local: via multi-turn potentiometer (kV ADJ) Remote: via external voltage source 0 to 10V (accuracy ± 1%)
Interlock:	Short to GND through a 12V lamp: HV/ON, OPEN:HV/OFF
Protection:	Over voltage, over current protection. Arc, short circuit.
Temperature Conditions:	Operational: 0 to 45°C Storage: -20 to + 85°C
Temperature Coefficient:	0.01 % per °C, voltage and current
Dimensions:	5.5" H x 3.3" W x 9.45" D (140mm x 83.5mm x 240mm)
Weight:	7.9 lbs. (3.6kg)
Regulatory & Safety:	Meets the requirements of IEC61010-1:2010, EN61010-1: 2010, UL61010-1: 2012, CA N/CSA C22.2 No. 61010- 1:2012 and 2006/95/EC Low Voltage Directive. Product carries the TUV SUD c/us mark.



Shasta 50kV & 60kV Power Supplies



Pin Assignments




Pin	Signal	Pin	Signal
J1	HV OUTPUT	J5	1 +10VDC
J2	1 +24VDC		2 N.C.
	2 GND		3 kV Prog Input
J3	1 FILAMENT		4 RMT/LCL kV Prog Selector
	2 FILAMENT RETURN (GND)		5 N.C.
	3 BIAS		6 mA Prog Input
J4	1 MONITOR RETURN (GND)		7 RMT/LCL mA Prog Selector
	2 kV MONITOR		8 N.C.
	3 mA MONITOR		9 GND
	4 INTLK		

Product Ordering Table

Part Number	Voltage	Power	Grid Control
9700001	50kV	50W	N/A
9700002	50kV	50W	0 to -300V
9700003	60kV	60W	N/A
9700004	60kV	60W	0 to -300V

Appendix A - Tubes



April 2016





Part Number	Outline Drawing	Target	Operating Range (kV)	Max Anode Current (mA)	Max Anode Power (W)	Max Filament Current (A)	Spot Size (μm)**
1000 Series Glass X-ray Tubes							
See data sheet on page 26 							
90501	8218	W	10 - 80	0.7	56	1.7	100 Max.
90502	8218	W	10 - 65	0.5	32.5	1.7	90 Max.
90503	8218	W	10 - 65	0.5	32.5	1.7	33 Nom.
90505	8225	W	10 - 50	1.0	50	1.7	200 Max.
90507	8218	W	10 - 80	0.5	40	1.7	33 Nom.
90512	8218	W	10 - 80	0.5	40	1.7	33 Nom.
1500 Series Glass X-ray Tubes with Window							
See data sheet on page 28 							
90011	8194	Rh	10 - 50	1.0	50	1.7	150 Max.
90030	8260	Mo	10 - 50	1.0	50	1.7	150 Max.
90034	8188	Cr	10 - 50	1.0	50	1.7	200 Typ.
90054	8205	W	10 - 50	1.0	50	1.7	150 Typ.
90066	8205	Mo	10 - 50	1.0	50	1.7	150 Typ.
90068	8188	Mo	10 - 50	1.0	50	1.7	150 Typ.
90069	8188	Cu	10 - 50	1.0	50	1.7	150 Typ.
90075	8214	Cu	10 - 30	1.0	50	1.7	150 Typ.
90077	8188	Rh	10 - 50	1.0	50	1.7	150 Typ.
90078	8205	Cr	10 - 50	1.0	50	1.7	200 Typ.
90083	8188	W	10 - 50	1.0	50	1.7	150 Typ.
90099	8189	W	10 - 50	1.0	50	1.7	70 Max.
90101	8205	Rh	10 - 50	1.0	50	1.7	150 Typ.
90112	8061	Mo	10 - 50	1.0	50	1.7	150 Typ.
90127	8268	Mo	10 - 50	1.5	75	1.7	150 Max.
90128	8268	W	10 - 50	1.5	75	1.7	150 Max.
90133	8260	Rh	10 - 50	1.0	50	1.7	150 Max.
90135	8188	Ag	10 - 50	0.5	25	1.7	150 Typ.
90140	8189	Cr	10 - 50	1.0	15	1.7	70 Max.
90169	8069	W	10 - 50	1.0	50	1.7	70 Max.
1501 Series Glass X-ray Tubes with Higher Current							
See data sheet on page 30 							
90015	8194	Rh	4 - 50	2.5	50	2.4	200 Max.
90122	8257	Cr	4 - 50	2.0	75	2.4	175 Max.
90131	8055	Rh	4 - 50	2.5	50	2.4	130 Max.
90132	8055	Mo	4 - 50	2.5	50	2.4	130 Max.

Part Number	Outline Drawing	Target	Operating Range (kV)	Max Anode Current (mA)	Max Anode Power (W)	Max Filament Current (A)	Spot Size (µm)**
1550 Series Glass X-ray Tubes with Grid Control							
See data sheet on page 32							
							
90200	8242	Cu	10 - 50	1.0	50	1.7	50 Max.
90201	8242	W	10 - 50	1.0	50	1.7	50 Max.
90202	8242	Mo	10 - 50	1.0	50	1.7	50 Max.
90203	8242	Co	10 - 50	1.0	10	1.7	50 Max.
90204	8242	Rh	10 - 50	1.0	50	1.7	50 Max.
90512	8218	W	10 - 80	0.7	56	1.7	33 Nom.
3000 Series Potted X-ray Tubes							
See data sheet on page 24							
							
90004	8145	W	4 - 30	0.5	15	1.3	1000 Typ.
90006	8053	W	4 - 30	0.2	6	1.3	1000 Typ.
90020	8052	Ti	4 - 30	0.5	15	1.3	1000 Typ.
90025	8052	W	4 - 10	1.0	10	1.3	1000 Typ.
90026	8263	Pd	4 - 30	0.1	2.5	1.3	1000 Typ.
90036	8053	W	4 - 30	0.3	9	1.3	1000 Typ.
90042	8195	Mo	4 - 30	0.5	15	1.3	1000 Typ.
90050	8053	W	4 - 30	0.2	6	1.3	1000 Typ.
90053	8204	Mo	4 - 15	1.0	15	1.3	1000 Typ.
90057	8156	Au	4 - 30	0.5	15	1.3	1000 Typ.
90072	8211	Cu	10 - 25	0.4	10	1.7	1000 Typ.
90116	8053	W	4 - 30	0.2	6	1.3	1000 Typ.
90118	8052	Fe	4 - 10	1.5	15	1.3	1000 Typ.
90125	8195	W	4 - 40	0.5	15	1.3	1000 Typ.
90136	8195	Cu	4 - 30	0.5	12	1.3	1000 Typ.
90141	8054	Rh	4 - 40	0.3	9	2.0	1000 Typ.
90142	8054	Ag	5 - 40	0.3	9	2.0	1000 Typ.
90143	8054	Mo	4 - 30	0.3	9	2.0	1000 Typ.
90144	8214	Cr	4 - 7.5	0.2	1.5	1.3	1000 Typ.
90145	8054	Pd	4 - 40	0.3	9	2.0	1000 Typ.
90146	8057	W	5 - 13.6	2.0	27.2	2.0	1000 Typ.
90151	8063	Rh	4 - 30	0.3	9	2.0	1000 Typ.
90152	8063	Ag	4 - 30	0.3	9	2.0	1000 Typ.
90153	8063	Mo	4 - 30	0.3	9	2.0	1000 Typ.

Appendix A - Tubes





April 2016

Part Number	Outline Drawing	Target	Operating Range (kV)	Max Anode Current (mA)	Max Anode Power (W)	Max Filament Current (A)	Spot Size (μm)**
Jupiter 5000 Series Radiation Shielded X-ray Tubes							
See data sheet on page 20 							
93000*	8166	W	10 - 50	1.0	50	1.7	165 Max.
93001	8166	Mo	10 - 50	1.0	50	1.7	150 Typ.
93009*	8166	Rh	10 - 50	1.0	50	1.7	150 Typ.
93048	8166	Cu	10 - 50	1.0	50	1.7	150 Typ.
93055	8166	Mo	10 - 50	1.0	50	1.7	150 Typ.
93057	8166	Rh	10 - 50	1.0	50	1.7	180 Typ.
93059	8203	Rh	10 - 50	1.0	50	1.7	180 Typ.
93069*	8166	W	10 - 50	1.0	50	1.7	70 Max.
93070	8166	Cr	10 - 50	1.0	50	1.7	200 Typ.
93071	8203	W	10 - 50	1.0	50	1.7	150 Typ.
93073	8166	Pd	10 - 50	1.0	50	1.7	200 Max.
93075*	8166	Rh	10 - 50	1.0	50	1.7	150 Max.
93078*	8203	Cu	10 - 50	1.0	50	1.7	175 Max.
93079*	8203	Mo	10 - 50	1.0	50	1.7	150 Typ.
93083*	8166	W	10 - 50	1.0	50	1.7	70 Max.
93087	8203	Rh	10 - 50	1.0	50	1.7	165 Max.
93089*	8166	W	10 - 50	1.0	50	1.7	50 Max.
93091*	8166	W	10 - 50	1.0	50	1.7	165 Max.
93093*	8167	W	10 - 50	1.0	50	1.7	70 Max.
93095*	8166	Mo	20 - 50	1.0	50	1.7	55 Max.
93111*	8232	W	10 - 50	1.5	75	1.7	150 Typ.
93115*	8232	Cu	10 - 50	2.5	50	2.5	150 Typ.
93119*	8252	Mo	10 - 50	1.5	75	1.7	150 Max.
93303*	8231	W	10 - 50	1.5	75	1.7	150 Typ.
93401	8166	W	10 - 50	1.0	50	1.7	150 Typ.
Jupiter 5000 Series Radiation Shielded X-ray Tubes with Brass Package							
See data sheet on page 20 							
93200	8235	Mo	10 - 50	1.0	50	1.7	150 Max.
93201	8235	W	10 - 50	1.0	50	1.7	175 Max.
93202	8235	Cr	10 - 50	1.0	50	1.7	175 Typ.
93203	8235	Mo	10 - 50	1.0	50	1.7	150 Max.
93208	8235	W	10 - 50	1.0	50	1.7	150 Max.

Part Number	Outline Drawing	Target	Operating Range (kV)	Max Anode Current (mA)	Max Anode Power (W)	Max Filament Current (A)	Spot Size (μm)**
Additional 5000 Series Radiation Shielded X-ray Tubes							
							
93117*	8232	Mo	4 - 50	1.5	75	1.7	170 Typ.
93025	8166	Ag	4 - 50	1.0	50	1.3	1000 Typ.
93006*	8232	Rh	4 - 50	2.5	50	2.5	150 Typ.
93033*	8203	Cu	4 - 50	1.0	50	2.4	250 Max.
93035	8166	Au	4 - 50	1.0	25	1.3	1000 Typ.
93046	8208	Mo	4 - 50	1.0	50	1.3	1000 Typ.
93066	8166	Rh	4 - 50	1.0	50	1.3	1000 Typ.
93072	8166	Ti	4 - 50	1.0	50	1.3	1000 Typ.
93080*	8232	Mo	4 - 50	2.5	50	2.5	150 Typ.
93084*	8166	Cr	4 - 50	2.5	50	2.4	150 Typ.
93138*	8270	Rh	4 - 50	2.0	50	2.5	1000 Typ.
93302*	8231	Mo	4 - 50	1.5	75	1.7	150 Typ.
93512*	8166	Fe	4 - 50	2.0	50	1.4	1000 Typ.
Neptune 5200 Series Water-Cooled Radiation Shielded							
See data sheet on page 18							
							
93211*	8250	Mo	10 - 50	2.0	100	2.4	175 Max.
93212*	8250	W	10 - 50	2.0	100	2.4	175 Max.
93215*	8250	W	10 - 50	2.0	100	2.4	175 Max.
93221*	8250	Rh	10 - 50	2.0	100	2.4	375 Max.
Apogee 5500 Series Radiation Shielded with Grid Control							
See data sheet on page 22							
							
93501*	8243	W	10 - 50	1.0	50	1.7	35 Nom.
93500*	8243	Cu	10 - 50	1.0	50	1.7	35 Nom.
93502*	8243	Mo	10 - 50	1.0	50	1.7	35 Nom.
93504*	8243	Rh	10 - 50	1.0	20	1.7	35 Nom.
93508*	8243	Cu	10 - 50	1.0	50	1.7	35 Nom.
93510*	8243	Cu	10 - 50	1.0	50	1.7	35 Nom.
93511*	8243	Mo	10 - 50	1.0	50	1.7	35 Nom.
Pinnacles 50 kV Microfocus Radiation Shielded X-ray Tube							
See data sheet on page 16							
							
Various - see data sheet	8062	W	10 - 50	1.0	12		10 Nom.

Appendix A - Tubes

April 2016

Part Number	Outline Drawing	Target	Operating Range (kV)	Max Anode Current (mA)	Max Anode Power (W)	Max Filament Current (A)	Spot Size (μm)**
Scafell Pike Integrated X-ray Source See data sheet on page 4 							
9100025	8100002	Rh	4 - 50	0.2	4		150 Max.
Trinity Integrated X-ray Source See data sheet on page 6 							
DS9100004	8300021	W	20 - 80	0.5	40		33 Max.
UltraBright Microfocus Source See data sheet on page 8 							
96002*	8236	Mo	20 - 60	2.0	60		20 Max.
96000*	8236	Cu	20 - 60	2.0	60		20 Max.
96004*	8236	W	20 - 90	2.0	80		20 Max.
96008*	8228	W	20 - 90	2.0	80		20 Max.
Nova Water-Cooled Microfocus Source See data sheet on page 12 							
96013*	8240	W	20 - 90	2.0	80		20 Max.
96016*	8240	Mo	20 - 60	2.0	60		20 Max.
Other Industrial X-ray Tubes							
90046	8162	W	4 - 50	5.0	250	2.4	500 Typ.
90098	8162	Fe	4 - 50	5.0	200	2.4	500 Typ.

* Includes a thermal switch which adds an additional level of protection to the cooling system safeguards.

** Inquire for outline drawings, detailed specifications and exact spot size dimensions.
 Max. = Maximum, Typ. = Typical, Nom. = Nominal (per IEC60336, NEMA XR5-1999)

Cable Part Number	Description	Length	Power Supply Part Number
Cables for Stainless Steel Packaged 5000 Series X-ray Tubes			
(All 5000 Series tubes unless noted below)			
			
9200008	Shasta High Voltage Cable	1m	Shasta 9700001
9200009	Shasta High Voltage Cable	2m	Shasta 9700001
9200010	Shasta High Voltage Cable	3m	Shasta 9700001
9200014	Shasta Low Voltage Cable	1m	Shasta 9700001
9200015	Shasta Low Voltage Cable	2m	Shasta 9700001
9200016	Shasta Low Voltage Cable	3m	Shasta 9700001
Cables for Apogee 5500 Series X-ray Tubes			
(All Apogee 5500 Series tubes)			
			
9200008	Shasta High Voltage Cable	1m	Shasta 9700002
9200009	Shasta High Voltage Cable	2m	Shasta 9700002
9200010	Shasta High Voltage Cable	3m	Shasta 9700002
9200011	Shasta Low Voltage Cable	1m	Shasta 9700002
9200012	Shasta Low Voltage Cable	2m	Shasta 9700002
9200013	Shasta Low Voltage Cable	3m	Shasta 9700002
Cables for Stainless Steel Packaged 5000 Series X-ray Tubes			
(Tube part numbers 93006, 93080, 93084, 93111, 93115, 93119, 93138, 93302, 93303, and 93512)			
			
92103	High Voltage Cable	1m	97013
92115	High Voltage Cable	2m	97013
92104	High Voltage Cable	3m	97013
9200014	Shasta Low Voltage Cable	1m	97013
9200015	Shasta Low Voltage Cable	2m	97013
9200016	Shasta Low Voltage Cable	3m	97013
Cables for Brass Packaged 5000 Series X-ray Tubes			
(Tube part numbers 93200, 93201, 93202, 93203, and 93208)			
			
92000	LGH High Voltage Cable A1	1m	SHASTA 9700001
92001	LGH High Voltage Cable A1	3m	SHASTA 9700001
9200014	Shasta Low Voltage Cable	1m	SHASTA 9700001
9200015	Shasta Low Voltage Cable	2m	SHASTA 9700001
9200016	Shasta Low Voltage Cable	3m	SHASTA 9700001
Cables for Neptune 5200 Series X-ray Tubes			
(All Neptune 5200 Series tubes)			
			
92103	LGH High Voltage Cable A1	1m	97013
92104	LGH High Voltage Cable A1	3m	97013
9200014	Shasta Low Voltage Cable	1m	97013
9200015	Shasta Low Voltage Cable	2m	97013
9200016	Shasta Low Voltage Cable	3m	97013

Application

For X-ray tubes that have been in storage or inactive for a period of three months or more.

Description

After a period of not being used, typically three months or more, residual gasses are released from the internal surfaces of the X-ray tube and accumulate into the tube vacuum. If the maximum rated voltage (in kV) is applied after a period of storage without performing a conditioning procedure, permanent damage to the X-ray tube may occur due to the destructive nature of high voltage arcs in the ionized gas. The following conditioning procedure is appropriate for both "new tubes", as there may have been a period of storage, as well as tubes that have been stored for three months or more.

Procedure

To prevent this damage the following special conditioning process should be followed:

- Adjust the kV to the lowest kV that your specific tube is rated for. Set the beam current to 0 mA and if any instability is noted on the mA meter allow it to stabilize to display 0mA. Operate at this condition for a minimum of 15 minutes.
- While maintaining the kV set in the previous step, adjust the beam current to 20% of rated maximum. Maintain this setting for 5 minutes or longer, until no instability is noted on the mA meter.
- Increase high voltage in 5kV steps at 5 minute intervals until 50% of maximum rated kV is reached. Hold 5 minutes at these conditions.
- Increase beam current to maximum rated mA.
- Continue to increase high voltage as before, in 5kV steps every 5 minutes until maximum rated kV or your maximum operating kV is reached. Allow at least 5 minutes at full power to insure that the tube is operating correctly in your system.

Note

If instability (especially loud popping) is observed, lower the kV setting to the previous step. Allow mA to stabilize for at least 5 minutes before increasing settings again.

Summary

Oxford Instruments offers X-ray tubes with different anode materials designed to suit a wide variety of applications. The anode material defines an X-ray tube’s characteristic spectrum. This application note shows the typical spectra of several different anode materials. The spectra provided are for reference only; your spectrum may differ from these according to the particular model of detector you are using, the geometry of your measurement setup, and the voltage and current on your X-ray tube.

X-ray Spectrum Theory


X-ray production involves bombarding a metal target in an X-ray tube with high-speed electrons that have been accelerated by tens to hundreds of kilovolts of electric potential. The electrons can eject other electrons from the inner shells of the atoms of the metal anode. Those vacancies will be filled when electrons drop down from higher energy levels and emit X-rays. These are known as characteristic X-rays and they have sharply defined energies associated with the difference between the atomic energy levels of the anode atoms. The Bohr atomic model predicts the energies of the characteristic X-rays. An X-ray spectrum is partially defined by the “peaks” or “lines” that result from bombarding different anode materials with highly accelerated electrons.

In addition to the characteristic peaks, an X-ray spectrum also has a background radiation pattern called the “Bremsstrahlung.” Bremsstrahlung means “braking radiation” and describes the radiation that is emitted when electrons are decelerated through a metal anode. The deceleration leaves behind excess energy, some of which is emitted in the form of radiation. Decelerated charges give off electromagnetic radiation, and when the energy of the electrons is high enough, that radiation is in the X-ray region of the electromagnetic spectrum.



Thus, the X-ray spectrum that is emitted from your X-ray tube is a combination of the characteristic peaks of the specific anode material and the Bremsstrahlung radiation that is present in all X-ray tubes.



Experimental Setup

The following spectra were gathered by pointing Oxford Instruments X-Ray Technology’s tubes directly at a Si-PIN photodiode detector system. There are a total of approximately one million counts in each spectrum.



Typical X-ray Spectra by Anode Material (cont.)





Summary

An X-ray tube is constrained in its operating range by four factors – maximum filament current, maximum power delivered to the anode, and maximum and minimum anode voltage. By operating your X-ray tube within these parameters, you may be able to achieve better results for your specific application while ensuring maximum longevity of your X-ray tube. This application note clarifies the constraints above and shows how an operating range is constructed. You can find all the particular values for your X-ray tube described in this application note on the datasheet.

Maximum Filament Current


The maximum filament current is a very strict constraint that prevents the filament from burning out, just like the filament in an incandescent light bulb. Like any other wire, a filament will melt because it cannot dissipate the heat generated from excessive current. Oxford Instruments X-Ray Technology has conducted extensive testing to determine the maximum current the filament in your X-ray tube can withstand. A common value is 1.7A, but this value varies by filament type and is given on the datasheet that comes with your X-ray tube. It constrains the first part of the operating range before the maximum power requirement takes over.

Maximum Power

Like the filament current limit, the power limit is a strict constraint that prevents the target from sublimating. An X-ray tube accelerates a very narrow beam of electrons to the target with a total power $P = IV$, where I is the beam current (not to be confused with the filament current – the current delivered to the filament itself) and V is the anode voltage. As you can see, this total power limit does not necessarily prevent using a higher beam current or voltage at a given power. Because beam current and voltage are inversely proportional in this relationship, raising one and lowering the other may still allow you to operate the X-ray tube within the maximum power constraint. Keep in mind that the anode voltage is limited as well, as detailed below. The maximum power constraint takes over after the filament limit is no longer a factor in the operating range.

Minimum and Maximum Anode Voltage

An X-ray tube requires a minimum high voltage applied to the anode in order to draw off electrons from the filament. When this condition is satisfied, the beam of electrons will form and accelerate towards the target. Below the minimum anode voltage, electrons will not be drawn off the filament, and thus the tube will produce no X-rays. At voltages lower than the minimum, some power supplies will overdrive (and potentially melt) the filament in an attempt to produce beam current when there are no electrons available. Our Shasta power supplies are designed to prevent damage to the filament. It is extremely important that you do not attempt to obtain beam current below this minimum anode voltage to avoid damaging the filament. On the other hand, the X-ray tube can only stand off a maximum high voltage applied to the anode. Beyond this voltage, arcing will occur and this can severely damage your X-ray tube. Both the minimum and maximum high voltages are sharp cut-offs that form the left and right edges of the operating range.



Conclusion

Some applications may require different settings than the typical “full power” at which most customers operate their X-ray tubes. By following the guidelines in this document, you may be able to achieve more desirable conditions for your application that still fall within the operating range constraints. In summary, to achieve the best possible conditions for your application, operate your X-ray tube within the constraints of maximum filament current, maximum power, and minimum and maximum anode voltage as described above.

When purchasing an X-ray tube, one of the most important questions which must to be answered is: how will the tube be packaged? System designers often put a lot of thought into the tube specifications, such as target material, spot size, etc, but the physical packaging of the tube can be a critical design choice. Many factors, including heat dissipation, radiation shielding, and design time must be considered.

Bare Tube

Bare tubes are just that – X-ray tubes with nothing else. It is incumbent on the system designer to design the radiation shielding, the insulating material, the high voltage and filament connections, and power supply integration. This can be quite a complex task, and is generally only appropriate for very large volume systems with specific requirements that cannot be met with Oxford's proven tube packaging solutions.

Potted Tube

Potted tubes are encapsulated in a silicone rubber material to provide electrical isolation and, in some cases, radiation shielding. High and low voltage cables may be included in the potting to aid connection to the X-ray power supply. Potted tubes provide an easier integration option than bare X-ray tubes however heat dissipation in a potted tube can be a challenge, and so potted tubes tend to be appropriate in low power or low duty cycle applications.

Packaged Tube

Packaged tubes are enclosed in a metal housing which acts as both a radiation shield and a cooling vessel. The packages are filled with a high dielectric liquid which both prevents high voltage breakdown (arcing) and effective cooling, requiring only an external fan to provide 50W of continuous power in many applications. Higher power packages with integrated water cooling systems are also available. Oxford's packaged tubes are fitted with connectors for easy plug-and-play operation with our Shasta X-ray power supply, enabling a quick setup procedure.

Integrated Source

Integrated X-ray sources include an X-ray tube, a high voltage and a low voltage power supply, and an analog or digital interface conveniently packaged in one box. This frees the system designer from all high voltage design concerns, and allows the X-ray device to be treated as a true "black box" component. Integrated solutions also speed up time-to-market, as the system designer only needs to integrate with a simple analog or digital interface, and won't be bogged down with often mysterious high voltage integration problems.

Many Oxford Instruments X-ray tubes come equipped with beryllium X-ray windows for maximum flux transmission. Beryllium is a metal that has low density and low atomic mass, and hence very low absorption of X-rays, making beryllium the preferred choice for X-ray tube windows where low energy transmission is desired.

Oxford Instruments also produces glass window tubes, which are much more robust than their beryllium counterparts, with the trade-off of decreased low energy flux. Glass window tubes are suitable for a wide variety of applications, including imaging and some types of analysis, and should be considered in harsh, humid, or debris filled environments.

If your X-ray tube has a beryllium window, please keep the following considerations in mind:

- The beryllium exit window is comprised of high purity vacuum tight beryllium metal, typically 127 microns thick.
- Beryllium can be toxic if improperly handled. Avoid contact with the beryllium window.
- The beryllium window is fragile and will be damaged by the slightest impact.
- Beryllium is highly soluble in polar solvents. Examples of polar solvents include water (including humidity), alcohol and acids. It is essential that you do not expose the beryllium window to these agents for prolonged periods of time, as they will destroy the beryllium window and compromise the internal high vacuum of the X-ray tube, causing it to fail.
- Unless absolutely necessary, all care should be taken to avoid any contact with the beryllium window, and tube installation should take into consideration keeping the window free of dust and debris. Should your beryllium exit window need to be cleaned, gently use a cotton swab and acetone (a non-polar solvent) and then immediately dry thoroughly with a cotton swab or soft dry air. Please note that damage to the beryllium window due to mishandling is not covered under your warranty.
- Helium is often used in X-ray spectroscopy. Helium is a very small atom and has a high transmission rate through the beryllium window. At a minimum, only beryllium exit windows of at least 127 microns should be considered when operating an X-ray tube in the presence of a Helium environment.
- If you operate an X-ray tube with a beryllium window in a vacuum environment, it is important to remember that the beryllium window is brittle and susceptible to damage caused by cycling between atmospheric pressure and vacuum environments typical for analytical analysis. Utilization of a secondary chamber is recommended to allow the X-ray tube to operate at sub-atmospheric pressures without cycling for each sample introduction.



Maximizing the Life of the X-ray Tube Filament

The process of producing electrons necessary for the production of X-rays in an X-ray tube begins by heating a tungsten wire. When heated to approximately 2000 degrees Celsius, tungsten is a copious emitter of electrons.

From this point several trade-offs in design become factors, which must be considered. The resulting design of a modern X-ray tube seeks to balance the relationship between performance and filament longevity.

Of importance to those users seeking a small X-ray focal spot, the relationship between a smaller wire filament and a small focal spot is well established. (This applies only to small focal spots when utilizing a tungsten wire filament. In the case of microfocus X-ray tubes, a dispenser cathode is typically employed.)

Since a smaller filament is preferred where possible, a typical filament “driver” circuit must be able to control the current to the filament quite carefully. This is due to the important relationship between filament current and actual temperature of the filament wire itself.

By example, the Jupiter Series 5000 X-ray tube requires more than 1.5 Amps current at 2 Volts to achieve the required filament temperature necessary for electron emission. However above 1.7 amps the filament enters a very high region of evaporation, and by 1.75 amps the filament reaches its melting point. Therefore careful control of the filament circuit is essential to a long lived X-ray tube. Our Shasta X-ray tube power supply has a tightly designed circuit which prevents the filament from exceeding its maximum allowable current. The Shasta power supply is perfectly matched to our X-ray tubes.

The process of heating the helical tungsten filament to produce electrons naturally causes the filament to evaporate. After a certain number of hours of normal operation, the filament will thin to the point of failure. The rate of filament evaporation, and thus the total number of hours required to thin the filament to the point of failure is a function of the chosen operating conditions.

The filament current required to heat and achieve a given X-ray beam current differs depending upon the required applied high voltage, as shown in Figure 1. To determine the anticipated life of a helical tungsten filament, one must estimate the average filament current employed throughout its life. Once estimated, the rate of evaporation can be used to estimate the normal filament life as shown in Figure 2. For example, if the user normally operates the X-ray tube at 40kV and 1.0 mA, this requires a filament current of approximately 1.60 A. Using the chart in Figure 2, this translates to approximately 40,000 hours of expected life.

A stand by condition of ~50% maximum filament current rating places the filament in a very low region of evaporation where the filament life is not measurably affected. You do not need to use a stand by condition to ensure maximum filament life, but you may find it beneficial as your power supply will achieve a steady state sooner.




Figure 1: Filament current required for the Jupiter Series 5000 X-ray tube




Figure 2: Filament life for the Jupiter Series 5000 X-ray tube

Summary

One of the most important safety considerations (along with the high voltage) in operating your X-ray tube from Oxford Instruments is how much shielding you should use to contain radiation. Shielding an X-ray tube involves a simple calculation based on mass attenuation coefficients for different materials, described below.

Disclaimer

Oxford Instruments does not make any claim that these calculated values will result in adequate attenuation. Due to material and geometry differences, these values may only be used as a starting point for your application. You must test your setup with a reliable dosimeter to ensure safety.

X-Ray Mass Attenuation Theory

This application note assumes a beam of monochromatic photons with an incident intensity I_0 that penetrates a material with mass thickness x and density. This beam will emerge with an intensity I given by the law¹

$$I/I_0 = \exp[-(\mu/x)x]$$

Values of μ/x have been empirically obtained using this equation and measured values of I_0 , I , and x . These values are compiled in the NIST X-Ray Mass Attenuation Coefficients² and are used for all the calculations in this document. Note that mass thickness is defined as the mass per unit area, and is obtained by multiplying the thickness t by the density such that $x = t$.

For composite materials such as Brass, the mass attenuation coefficients are obtained using a weighted average:

$$\mu/x = \sum w_n(\mu/x)_n$$

Here w_n is the fraction by weight of the n^{th} element in the material and, similarly, $(\mu/x)_n$ is the mass attenuation coefficient of the n^{th} element in the material.


We have completed the shielding calculations with a simple spreadsheet application for various materials using the theory above. The table below shows the shielding values for various materials. (Note that the Brass in this calculation is composed of 65% Copper, 33.5% Zinc, and 1.5% Lead). As a secondary consideration, these values have been crosschecked using existing Oxford Instruments XT experimental equipment.

Material	50 keV	100 keV
Lead	1.5 mm	3.0 mm
Iron	5.0 mm	20.0 mm
Brass	3.0 mm	12.0 mm


The following charts show the transmission characteristics by varying material thickness.




Transmission for Various Lead Thicknesses



Transmission for Various Iron Thicknesses





Conclusion

To ensure safety, it is extremely important to adequately shield the outside environment from X-rays that are being emitted from your X-ray tube. In order to do this, we recommend starting with at least the thickness of materials shown above and measuring the output with a dosimeter. Again, while these values have been developed from first principles, it is imperative that you measure any setup thoroughly before putting it into full use.

References

- ¹ <http://www.physics.nist.gov/PhysRefData/XrayMassCoef/chap2.html>
- ² <http://www.physics.nist.gov/PhysRefData/XrayMassCoef/tab3.html>

The most frequent mode of failure of X-ray tubes is the failure to adequately dissipate the heat generated during normal operation.

Greater than 99% of the kinetic energy imparted on the electron beam is lost in the form of heat at the anode target. Thus, a 50W X-ray tube will produce roughly 49.8W of energy in the form of heat just through the conversion process. Add to this the thermal energy produced by the helical tungsten filament and one can readily see that heat dissipation is a major factor.

Inadequate cooling of an X-ray tube can cause it to fail in two ways.

The first is sublimation of the anode target material. In converting the anode target material directly from a solid to a gas (sublimation), the resulting vapor rapidly degrades the internal high vacuum necessary for proper operation on the X-ray tube. The loss of high vacuum results in a failure of the X-ray tube to withstand the high voltage gap between the cathode electron source (helical tungsten filament) and the target anode. The X-ray tube begins to short circuit, or arc, which in turn liberates more gas that further degrades the internal vacuum, resulting in an X-ray tube that no longer functions.

The second failure mode caused by improper heat dissipation is the liberation of damaging ions. If the X-ray tube anode is allowed to surpass the vapor pressure point of the target material, ions will be liberated. These ions are attracted back toward the helical tungsten filament and begin to erode the filament through an ion scrubbing process. This can cause the filament to break, creating an open circuit.

Prevention of both of these failure modes is made possible by ensuring that the X-ray tube is not allowed to overheat. This means careful monitoring of the cooling system with fault protection in the event of a cooling system failure.

Many of our packaged X-ray tubes offer an integrated thermal switch that helps prevent permanent damage to the X-ray tube.

Please confirm that your cooling system can maintain the temperature range that is recommended on our products.




Image credits

1. Image courtesy of Molecubes (*see catalog back cover*)

X-ray Fluorescence Periodic Table

1 1.01 Hydrogen H	2 4.00 Helium He	3 6.94 Lithium K α : 0.05 L α : 0.26 Li	4 9.01 Beryllium K α : 0.11 L α : 0.34 Be	5 10.81 Boron K α : 0.18 L α : 0.28 B	6 12.01 Carbon K α : 0.28 L α : 0.28 C	7 0.39 Nitrogen K α : 0.39 L α : 0.39 N	8 15.99 Oxygen K α : 0.52 L α : 0.52 O	9 18.99 Fluorine K α : 0.68 L α : 0.68 F	10 20.18 Neon K α : 0.85 L α : 0.85 Ne	11 22.99 Sodium K α : 3.31 L α : 0.26 Na	12 24.31 Magnesium K α : 1.25 L α : 0.34 Mg	13 26.98 Aluminum K α : 1.49 L α : 0.28 Al	14 28.09 Silicon K α : 1.74 L α : 0.28 Si	15 30.97 Phosphorus K α : 2.01 L α : 0.28 P	16 32.06 Sulfur K α : 2.31 L α : 0.28 S	17 35.45 Chlorine K α : 2.62 L α : 0.28 Cl	18 39.95 Argon K α : 2.96 L α : 0.28 Ar	19 39.10 Potassium K α : 3.31 L α : 0.26 K	20 40.08 Calcium K α : 3.69 L α : 0.26 Ca	21 44.96 Scandium K α : 4.09 L α : 0.40 Sc	22 47.88 Titanium K α : 4.51 L α : 0.45 Ti	23 50.94 Vanadium K α : 4.95 L α : 0.51 V	24 51.99 Chromium K α : 5.41 L α : 0.57 Cr	25 54.94 Manganese K α : 5.89 L α : 0.64 Mn	26 55.85 Iron K α : 6.40 L α : 0.70 Fe	27 58.93 Cobalt K α : 6.93 L α : 0.78 Co	28 58.71 Nickel K α : 7.47 L α : 0.85 Ni	29 63.55 Copper K α : 8.04 L α : 0.93 Cu	30 65.37 Zinc K α : 8.63 L α : 1.01 Zn	31 69.72 Gallium K α : 9.24 L α : 1.10 Ga	32 72.59 Germanium K α : 9.88 L α : 1.19 Ge	33 74.92 Arsenic K α : 10.53 L α : 1.28 As	34 78.96 Selenium K α : 11.21 L α : 1.38 Se	35 79.90 Bromine K α : 11.91 L α : 1.48 Br	36 83.80 Krypton K α : 12.63 L α : 1.59 Kr	37 85.47 Rubidium K α : 13.28 L α : 1.69 Rb	38 87.62 Strontium K α : 14.14 L α : 1.81 Sr	39 88.91 Yttrium K α : 14.97 L α : 1.92 Y	40 91.22 Zirconium K α : 15.75 L α : 2.04 Zr	41 92.91 Niobium K α : 16.58 L α : 2.17 Nb	42 95.94 Molybdenum K α : 17.45 L α : 2.29 Mo	43 98.91 Technetium K α : 18.33 L α : 2.42 Tc	44 101.07 Ruthenium K α : 19.24 L α : 2.54 Ru	45 102.91 Rhodium K α : 20.12 L α : 2.67 Rh	46 106.4 Palladium K α : 21.18 L α : 2.80 Pd	47 107.87 Silver K α : 22.16 L α : 2.93 Ag	48 112.40 Cadmium K α : 23.17 L α : 3.06 Cd	49 114.82 Indium K α : 24.27 L α : 3.19 In	50 118.69 Tin K α : 25.27 L α : 3.32 Sn	51 121.75 Antimony K α : 26.36 L α : 3.45 Sb	52 127.60 Tellurium K α : 27.47 L α : 3.57 Te	53 126.90 Iodine K α : 28.46 L α : 3.69 I	54 131.30 Xenon K α : 29.54 L α : 3.81 Xe	55 132.91 Caesium K α : 30.65 L α : 3.94 Cs	56 137.34 Barium K α : 31.77 L α : 4.07 Ba	57 138.91 Lanthanum K α : 32.91 L α : 4.20 La	58 140.12 Cerium K α : 34.06 L α : 4.33 Ce	59 140.91 Praseodymium K α : 35.20 L α : 4.46 Pr	60 144.24 Neodymium K α : 36.36 L α : 4.59 Nd	61 144.91 Promethium K α : 37.51 L α : 4.72 Pm	62 150.36 Samarium K α : 38.73 L α : 4.92 Sm	63 151.96 Europium K α : 39.90 L α : 5.05 Eu	64 157.25 Gadolinium K α : 41.08 L α : 5.18 Gd	65 158.93 Terbium K α : 42.27 L α : 5.31 Tb	66 162.50 Dysprosium K α : 43.48 L α : 5.44 Dy	67 164.93 Holmium K α : 44.69 L α : 5.57 Ho	68 167.26 Erbium K α : 45.91 L α : 5.70 Er	69 168.93 Thulium K α : 47.14 L α : 5.83 Tm	70 173.04 Ytterbium K α : 48.37 L α : 5.96 Yb	71 174.97 Lutetium K α : 49.61 L α : 6.09 Lu	72 175.04 Hafnium K α : 50.91 L α : 6.19 Hf	73 178.49 Tantalum K α : 52.20 L α : 6.32 Ta	74 183.85 Tungsten K α : 53.49 L α : 6.45 W	75 186.21 Rhenium K α : 54.79 L α : 6.58 Re	76 193.22 Osmium K α : 56.08 L α : 6.71 Os	77 193.22 Iridium K α : 57.38 L α : 6.84 Ir	78 195.09 Platinum K α : 58.68 L α : 6.97 Pt	79 196.97 Gold K α : 59.98 L α : 7.10 Au	80 200.59 Mercury K α : 61.27 L α : 7.23 Hg	81 204.37 Thallium K α : 62.57 L α : 7.36 Tl	82 207.19 Lead K α : 63.87 L α : 7.49 Pb	83 208.98 Bismuth K α : 65.17 L α : 7.62 Bi	84 210 Polonium K α : 66.47 L α : 7.75 Po	85 210 Astatine K α : 67.77 L α : 7.88 At	86 222 Radon K α : 69.07 L α : 8.01 Rn	87 223 Francium K α : 70.37 L α : 8.14 Fr	88 226 Radium K α : 71.67 L α : 8.27 Ra	89 227 Actinium K α : 72.97 L α : 8.40 Ac	90 232.04 Thorium K α : 73.87 L α : 8.53 Th	91 231 Protactinium K α : 74.77 L α : 8.66 Pa	92 238.03 Uranium K α : 75.77 L α : 8.79 U	93 237 Neptunium K α : 76.67 L α : 8.92 Np	94 244 Plutonium K α : 77.67 L α : 9.05 Pu	95 243 Americium K α : 78.57 L α : 9.18 Am	96 247 Curium K α : 79.47 L α : 9.31 Cm	97 247 Berkelium K α : 80.37 L α : 9.44 Bk	98 251 Californium K α : 81.27 L α : 9.57 Cf	99 254 Einsteinium K α : 82.17 L α : 9.70 Es	100 257 Fermium K α : 83.07 L α : 9.83 Fm	101 256 Mendelevium K α : 83.97 L α : 9.96 Md	102 255 Nobelium K α : 84.87 L α : 10.09 No	103 257 Lawrencium K α : 85.77 L α : 10.22 Lr	104 258 Ununquadium K α : 86.67 L α : 10.35 Uuq	105 258 Ununpentium K α : 87.57 L α : 10.48 Uup	106 259 Ununhexium K α : 88.47 L α : 10.61 Uuh	107 259 Ununseptium K α : 89.37 L α : 10.74 Uus	108 262 Ununoctium K α : 90.27 L α : 10.87 Uuo
-----------------------------------	----------------------------------	---	---	---	--	---	--	--	--	--	---	--	---	---	---	--	---	--	---	--	--	---	--	---	--	--	--	--	--	---	---	--	---	--	--	---	--	---	--	--	---	---	---	---	--	--	---	--	---	--	---	---	---	---	--	---	--	--	---	--	--	--	--	---	--	---	--	---	---	--	---	--	---	---	--	---	--	--	---	--	--	---	---	---	--	---	---	---	---	---	--	--	--	--	---	--	--	--	---	---	---	---	---	---	--	---	--



X-RAY

Oxford Instruments X-Ray Technology

X-ray Fluorescence



Metal Analysis and Identification

Industrial Inspection



Printed Circuit Board Drilling and Inspection

Medical Imaging



Bone and Tissue Imaging



Energy and Mining



Food and Packaging Safety



Small Animal and Drug Studies¹

Visit www.oxford-instruments.com/xt or xray-sales@oxinst.com for more information

This publication is the copyright of Oxford Instruments plc and provides outline information only, which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or regarded as the representation relating to the products or services concerned. Oxford Instruments' policy is one of continued improvement. The company reserves the right to alter, without notice the specification, design or conditions of supply of any product or service. Oxford Instruments acknowledges all trademarks and registrations. © Oxford Instruments plc, 2016. All rights reserved. Document reference: Part no: 3900003 - April 4, 2016



The Business of Science®



X-Ray Technology
360 El Pueblo Road, Suite 104
Scotts Valley, CA 95066, USA

Phone: +1 (831) 439-9729
Fax: +1 (831) 439-6050
Email: xray-sales@oxinst.com